DECUS PERSONAL

Stimulator za celjenje ran

OPIS - NAVODILA ZA UPORABO STIMULATORJA ZA CELJENJE RAN
--

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Literatura

1. Dinsdale S.M.: Decubitus ulcers: Role of pressure and friction in causation, *Arch. Phys. Med.*, 55:147-152, 1974.
2. Brand P.W.: Pressure sore: The problem, *Bedsore biomechanics*, pp. 19-25, 1975.
3. Krouskop T.A.: A synthesis of the factors that contribute to pressure sore formation, *Medical Hypotheses*, 11:255-267, 1983.
4. Mathews A.P.: Electrical polarity in the hydroids, *Amer. J. Physiology*, 8:294-299, 1903.
5. Assimacopoulos D.: Low intensity negative electric current in the treatment of ulcers of the leg due to chronic venous insufficiency, *Am.J. of Surg.*, 115:683-687, 1968.
6. Wheeler P.C., Wolcott L.E., Morris J.L., Spangler M.R.: Neural considerations in the healing of ulcerated tissue by clinical electrotherapeutic application of weak direct current: Findings and theory. *Proc. Neuroelectric Conf.* (Ed. CC Thomas), pp: 83-99, 1970.
7. Edel H., Freund R.: Gleichstrombehandlung chronischer Hautulzera und sekundär heilender Wunden, *Zschr. Physiotherapie*, 27:457-464, 1975.
8. Thomas D., Lassaux G., Bacquillon S., Daumard N.: Traitement d'une escarre du decubitus par electrostimulation, Premier resultat a propos d'un cas, *Cah. Kinesiother.*, 93:57-63, 1982.
9. Carley P.J., Wainapel S.F.: Electrotherapy for acceleration of wound healing: Low intensity direct current, *Arch. Phys. Med Rehabil.*, 66:443-446, 1985.
10. Alon G.: Healing currents?, *IEEE Spectrum*, p. 22, May 1987.
11. Tallis R.C., Illis L.S., Sedgwick E.M., Hardwdige C., Kennedy K.: the effect of spinal cord stimulation upon peripheral blood flow in patients with chronic neurological disease, *Int. Rehab. Med.*, 5:4-9, 1983.
12. Levine S.P., Caderna P.S., Brooks S.V., Freedman R.H.: Functional electrical stimulation for the prevention of pressure sores, *Proc. IEEE Seventh Ann. Conf. Engn in Med and Biol. Soc.*, 707-710, 1985.

13. Maležič M., Benko H., Vodovnik L., Stefanovska A.: Dekubitus - Celjenje ran z električno stimulacijo, Zbornik del XXX Jugosl. Konf. ETAN, Herceg-novi, XI.93-100, 1986.
14. Benko H., Maležič M., Košorok V., Vodovnik L., Stefanovska A., Turk R., Dimnik Z., Obreza P.: Our experiences with healing of wounds with the use of electrical stimulation, Seminar of functional electrical stimulation, ZRI/ Ljubljana, pp.94-100, October 1986.
15. Vodovnik L., Stefanovska A., Benko H., Maležič M., Turk R., Košorok V.: Soft Tissue healing with low frequency pulsed currents, Trans. Sixt Ann. Meeting BRAGS, Utrecht, p.74, 1986.
16. Stefanovska A., Vodovnik L., Benko H., Maležič M., Turk R., Košorok V.: Regeneration of ulcerated tissue by electrical stimulation, Proc. 2nd Vienna Int. Workshop on Funct. Electrostimulation, Vienna, pp. 307-310, 1986.
17. Stefanovska A., Vodovnik L., Benko H., Maležič M., Turk R., Kolenc A., Turk R.: Enhancement of ulcerated tissue healing by electrical stimulation, Proc. 10th Annual Conf. RESNA, pp.585- 587, 1987.
18. Stefanovska A., Vodovnik L., Benko H., Kolenc A., Maležič M., Turk R.: Exogenous electrical currents: Influence on wound healing, Trans Seventh Ann. Meeting BRAGS, p.22, 1987.
19. Košorok V., Benko H., Vodovnik L.: Zdravljenje preležanin s pomočjo električne stimulacije, Zbornik del 5 Kongresa lekara za fizikalnu medicinu i rehabilitaciju Jugoslavije, st. 35-38, 1987.

Institut "Jožef Stefan", Ljubljana, Slovenija

Jamova c. 39, 1000 Ljubljana

Telefon: 01 477 3400 e-mail: stimulatorji@ijs.si