

MICROFES – Electrical stimulator for drop-foot correction

TECHNICAL DATA		
operating mode	automatic adjustment to pace	
pulses	biphasic charge balanced current stimulating pulses	
amplitude	mA	0 – 50
frequency	Hz	25
pulse width	ms	0,25
stimulation-on time	s	0 -5
stimulation-off delay after heel contact	ms	50 – 350, adjusted by the clinican
power consumptin	mA	30 (max)
size	mm	64 x 46 x 20
weight (including battery)	g	140
power source	1 x 1,5 V IEC LR6 (alkaline cell)	