

ALT-2 – Alternating two channel electrical stimulator

TECHNICAL DATA		
operating mode	Alternating (4s/4s) intermittent trains of stimulating pulses, or external manual switch control, two channels	
pulses	biphasic charge balanced current stimulating pulses	
amplitude	mA	0 – 180
frequency	Hz	20
pulse width	ms	0,25
power consumption	mA	200 (max)
size	mm	100 x 57 x 24
weight (including batteries)	g	120
power source	2 x 1,5 V IEC LR6 (alkaline cell) 2 x 1,2 V IEC HR6 (NiMH rechargeable cell)	