

NEVROBOL Personal – Electrical stimulator for pain relief

TECHNICAL DATA		
operating mode	continious train of stimulating pulses, or external manual switch control	
pulses	biphasic charge balanced current stimulating pulses	
amplitude	mA	0 – 50
frequency	Hz	100
pulse width	ms	0,3
power consumption	mA	55 (max)
size	mm	120 x 65 x 30
weight (including batteries)	g	140
power source	2 x 1,5 V IEC LR6 (alkaline cell) 2 x 1,2 V IEC HR6 (NiMH rechargeable cell)	